THERMO SLEEVE

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THERMO SLEEVE

PRODUCT DESCRIPTION	22
INTENDED USE	23-24
SAFETY INSTRUCTIONS	25
GENERAL APPLICATION TIPS FOR USE	26
WHY USE HEAT?	27-28
INSTRUCTIONS FOR HEAT THERAPY	29
WHY USE COLD?	30-31
INSTRUCTIONS FOR COLD THERAPY	32
INSTRUCTIONS FOR APPLYING THE SHOULDER SLEEVE	33
GENERAL TIPS	34
SIZING TABLE	35-36
SYMBOLS	37
WARRANTY DECLARATION	38

CONTENT

PRODUCT DESCRIPTION

The Reboots Thermo Sleeves are reusable products for localized heat and cold applications in the area of muscles, joints and limbs in the sports, fitness and recovery sector.

The high-quality gel pad in soft polyester enables even 360° heat or cold application, adapts flexibly and sits comfortably on the affected area - without dripping or slipping.

The sleeves are suitable for use as part of your personal regeneration routine, for example after exercise or when you need precise heating or cooling.

The application is exclusively for your well-being and does not replace medical treatment.

THE SLEEVES CONSIST OF:

- a flexible gel pad for temperature conduction
- a soft and flexible, textile outer cover
- (depending on the model) an elastic or adjustable fixing element

THE SLEEVES ARE AVAILABLE IN THE FOLLOWING VERSIONS:

- Thermo Sleeve (arm/leg) - for elbows and knees
- Thermo Sleeve Full Leg - from the ankle to the groin
- Thermo Sleeve Ankle - for ankle and foot area
- Thermo Sleeve Shoulder - for shoulder and upper back/chest area

INTENDED USE

The Reboots Thermo Sleeves are designed for the external application of heat or cold to specific areas of the body.

They are intended to support the subjective feeling of well-being after sporting or physical exertion.

The products are not intended for the diagnosis, treatment, cure or prevention of disease.

Reboots Thermo Sleeves are not a medical product. They are not a substitute for medical treatment. Medical advice should be sought if symptoms persist.

USE AS DIRECTED:

- External use only
- Use on intact skin
- Use according to the temperature specifications
- Adhere to the maximum recommended duration of use
- Use by adults or under supervision

NOTE FOR USE IN THE SHOULDER AREA:

- Do not fall asleep lying on the sleeve
- Do not use if you have an existing cardiovascular disease without consulting a doctor
- Do not place over large areas of the chest and collarbone

TARGET GROUP:

- Athletically active persons
- People with a need for localized heat or cold application as part of personal regeneration

NOT SUITABLE FOR:

- Children under 12 years of age
- Persons with severe circulatory disorders
- People with impaired temperature perception
- Persons with diabetes, cardiovascular diseases, impaired temperature perception and pregnant women should seek medical advice before use.

FORESEEABLE MISUSE:

The following in particular are considered improper use:

- Use while you are asleep
- Use with reduced temperature sensitivity
- Application to open wounds
- Heating for longer than the recommended time
- Fixing too tightly so that breathing is restricted (shoulder)
- Application directly over implanted electronic devices
- Permanent cooling below 0°C directly on the skin

SAFETY

INSTRUCTIONS

DO NOT USE FOR:

- known cardiovascular diseases without consulting a doctor
- pronounced circulatory disorders
- Sensory disorders in the shoulder/arm area
- severe neurological diseases
- acute inflammation with fever
- open wounds
- Implanted electronic devices in the chest area (e.g. pacemaker)

GENERAL APPLICATION TIPS

TO ENSURE THAT YOU ENJOY THE BEST POSSIBLE EXPERIENCE WITH THE THERMO SLEEVES, PLEASE NOTE THE FOLLOWING INSTRUCTIONS:

- **Practical storage & cooling:** The additionally available Thermobag serves as a practical transport and cooling bag so that your sleeve is always ready for use.
- **Dry use:** Make sure that the Thermo Sleeve does not get wet if possible to preserve the material quality in the long term.
- **Individual sensitivity to cold:** Everyone feels cold differently. Adjust the cooling time individually so that it is comfortable for you.
- **Safety first:** If you are unsure how to use the Thermo Sleeve in your particular case, consult your doctor to ensure the best use for your needs.

WHY USE HEAT?

WHY USE HEAT AFTER TRAINING?

Many athletes find heat pleasant and relaxing after physical exertion. The application of heat can help you to calm down after intensive sessions and positively support your own body awareness.

WHY USE WARM COMPRESSES?

Applying heat to painful or tense areas can promote blood circulation, loosen muscles and thus provide deep comfort:

- Relieve muscle tension - Gently relaxes stiff, overworked muscles.
- Improve circulation - Promotes better blood flow for faster recovery.
- Relieves cramps and pain - Helps with joint pain and everyday discomfort.
- Promotes mobility - Facilitates movement by loosening tense areas.
- Promotes relaxation - A calming, cozy way to de-stress and unwind.

APPLICATION AREAS	EFFECTS	SAFETY INSTRUCTIONS
after intensive training sessions	Pleasant sensation of warmth	Do not apply to open or irritated skin
for subjective feelings of muscle stiffness	can contribute to relaxation after training	
For relaxation after physical exertion	is often perceived as beneficial after sporting activity	Do not use if your temperature is restricted
		Discontinue use if it is perceived as unpleasant Discomfort/ shivering
		If unsure, seek medical advice before use

INSTRUCTIONS FOR HEAT THERAPY

MICROWAVE

- 1. Take the thermo sleeves out of the packaging to warm them up.
- 2. Place the sleeves in the microwave for about 20 seconds
- 3. **NOTE:** Depending on which microwave you use, it may take longer or shorter to warm up. It is best to heat the sleeve in 10-second increments.
- 4. Check that the sleeve is at the right temperature for you. It should feel pleasantly warm, not hot.
- 5. Repeat this step until the sleeve is at the right temperature for you.
- 6. **Important:** Do not heat your Thermo Sleeves for longer than 1 minute at a time and no longer than 2 minutes in total in the microwave!
- 7. If your Thermo Sleeves are frozen, allow them to thaw at room temperature before heating.

PERFORMANCE	TIME		
	S/M	L	XL
650-799 W	max. 60 seconds	max. 80 seconds	max. 120 seconds
800-1000 W	max. 40 seconds	max. 60 seconds	max. 90 seconds

Only heat the Thermo Sleeves from room temperature and for optimum results, it is best to heat them in 10-second increments.

WHY

USE COLD?

WHY USE COLD AFTER TRAINING?

In the sports and fitness sector, cold is often perceived as refreshing and invigorating - especially after intensive exertion or long training sessions. Many users use cold applications to achieve a pleasant feeling of freshness after training.

WHY COLD COMPRESSES?

A cold application has an anti-inflammatory, pain-relieving and circulation-promoting effect and thus contributes to this:

- Minimize swelling - Cold constricts blood vessels, reducing excess fluid build-up.
- Alleviates discomfort and sensitivity - Helps to numb discomfort and relieve sensitivity.
- Calms inflammation - Reduces inflammatory reactions that contribute to flare-ups.
- Promotes circulation and lymph flow - Stimulates the removal of fluids.
- Instant relief - A quick and easy way to feel better whenever you need it.

Integrating cold therapy into your routine can make a big difference to your daily wellbeing and mobility.

APPLICATION AREAS	EFFECTS	SAFETY INSTRUCTIONS
after physical exertion	Cooling, refreshing feeling	Do not apply to open or irritated skin
For a cooling feeling of freshness	supports the personal regeneration routine	Do not use if your temperature is restricted
as part of your personal recovery routine		Discontinue use if it is perceived as unpleasant Discomfort/ shivering
		If unsure, seek medical advice before use

REFRIGERATOR

1. The Thermo Sleeve should not be too cold for cold applications. A temperature of 2-8 degrees is ideal.
2. Place your sleeve in a resealable, suitable bag in the refrigerator.
3. Chill the sleeve for at least two hours.
4. You can simply store the thermo sleeve in the refrigerator until you are ready to use it.
5. Do not use the cooled sleeve in its intended place for longer than 15 minutes.

NOTE: If you want to cool down immediately after your workout on the go, simply take your cooled Thermo Sleeve with you in a small cool bag.

FREEZING CABINET

1. Need to be quick? Then you can also cool your Thermo Sleeves in the freezer.
2. Place your sleeve in a resealable bag (e.g. in the bag provided or another sealable container) in the freezer.
3. Cool the sleeve for a maximum of 20 minutes.

NOTE: Make sure that the sleeve does not get too cold in the freezer. The cooling should still feel pleasant.

NOTE: Take a break of at least 30 minutes between two applications.


INSTRUCTIONS FOR CREATING

1. Place the sleeve over the affected shoulder
2. Position the gel pad centered on the shoulder area.
3. Guide the belt from the back to the front.
4. Secure so that the sleeve fits snugly but does not constrict.


Note: Breathing must not be restricted

1

2

3

4

GENERAL TIPS

- 1** We recommend using a sealable bag (e.g. our separately available Thermobag) to store the sleeve.
- 2** Try not to get your thermal sleeves wet.
- 3** Everyone feels differently. Even if we give you recommendations for use here, please listen to your own feelings of warmth and cold and warm/cool the sleeve so that it is comfortable for you.
- 4** If you are unsure about the best way to use the Thermo Sleeve in your case, it is always best to consult your doctor.

SIZE
TABLE

SHORT

	S	M	L
Product dimensions	26cm length 23cm circumference	26cm length 33cm circumference	26cm length 53cm circumference
Appropriate scope	26-35cm	33-44cm	50-63cm

FULL LEG

	ONE SIZE
Product dimensions	Length: 75 cm; top circumference: 34 cm bottom circumference: 16 cm
Appropriate scope	30-44cm at the top of the thigh

ANKLE


	ONE SIZE
Product dimensions	35x20cm
Appropriate scope	Ankle circumference: 17.78-25.4cm
Weight	580g

SHOULDER

	ONE SIZE
Product dimensions	35x20cm
Appropriate scope	Chest circumference: 74 - 95cm Biceps circumference: 29 - 40cm
Weight	Weight1000g

PRODUCT INFORMATION

Material: 85% nylon, 15% spandex

	Manufacturer
	Country of manufacture: China
	France: Recycle packaging accordingly
	Importer
	Follow the instructions for use
	Do not wash

	Do not bleach
	Hand wash only
	Do not tumble dry
	Non-ironing
	Dispose of the product in accordance with local regulations.
	RoHS compliant

SYMBOLS

WARRANTY DECLARATION

We grant 24 months on the product.

Warranty conditions: Insofar as there is a legal obligation to provide a warranty, either a replacement will be provided free of charge or the manufacturer will be given the opportunity to rectify the defect, excluding the right to rescission or reduction. If the repair is unsatisfactory despite several attempts or if a replacement device also has a defect for which the manufacturer is responsible, the customer has the right to cancel the contract or reduce the purchase price.

Exclusion of warranty: Excluded from the warranty are defects caused by violent damage, improper operation, external force or by modifications and repairs carried out by third parties, such as defects caused by incorrectly dimensioned or short-circuited fuses or defects caused by normal wear and tear.

Warranty requirement: The warranty claim is only valid in conjunction with the original proof of purchase (invoice). Therefore, keep the proof of purchase in a safe place. If you have any problems or questions when using the system, please do not hesitate to contact us on **+49 (0) 157 359 900 81** or **hallo@reboots.de**.

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Wellcosan GmbH, c/o Caya Postbox 794627,